

19:00


Traditional Kinetic Limulus Amebocyte Lysate (LAL) Assay Procedure Quick Guide

This is a step by step guide depicting how to perform a traditional kinetic LAL assay. Prior to initiating the assay procedure, allow reagent vials to equilibrate to room temperature. The incubating microplate reader should also be turned on and a plate template created in the WinKQCL® Software.

Step 3


Label the tubes with the appropriate endotoxin concentration and add 0.9 mL of LRW to each.

(Example based on a test with an operating standard curve of 0.005-50 EU/mL.)




Step 5

Dispense 100 µL of the LRW blank, endotoxin standards, product samples, positive controls, etc. into the appropriate wells of the microplate.




Step 8

Pour LAL into a reagent reservoir and mix gently.




Step 1

Reconstitute Control Standard Endotoxin (CSE) with LAL Reagent Water (LRW) to yield a solution containing 50 EU/mL or 100 EU/mL depending on assay method being used.




Step 4

Prepare a series of endotoxin standards.



Step 6

Pre-incubate the plate for ≥10 minutes at 37°C ± 1°C in the microplate reader.



Step 9


Use an eight channel pipettor to dispense 100 µL of LAL into the appropriate wells of the microplate.



Step 2

Vortex for 15 minutes.

Immediately prior to use, reconstitute LAL and gently swirl.



Step 10

Initiate the test by clicking the OK button in the WinKQCL® Software.



Materials, equipment and documents needed

Reagents

- Limulus Amebocyte Lysate (LAL) Reagent (Kinetic-QCL® or PYROGENT® 5000 Reagent)
- Control Standard Endotoxin (CSE)
- LAL Reconstitution Buffer (Required for the PYROGENT® 5000 Kinetic Turbidimetric LAL Assay)
- LAL Reagent Water (LRW) (# W50-640, W50-100, W50-500)

Kits are available in a wide range of sizes and configurations. Please contact your local sales representative for additional information.

Accessories

- Glass dilution tubes (# N207)
- Individually wrapped serological pipettes (optional)
- Tips
- 96-well plates (# 25-340)
- Reagent reservoirs (# 00190035)

Use pyrogen-free accessories that have been qualified for endotoxin testing.

Equipment and software

- Eight channel and pipettor
- Incubating absorbance microplate reader
- WinKQCL® Software
- Pipettors
- Timer
- Vortex mixer

Supporting documents

- Certificate of Analysis (CoA), www.lonza.com/coa
- Limulus Amebocyte Lysate (LAL) Kinetic-QCL® Package Insert or Limulus Amebocyte Lysate (LAL) PYROGENT® 5000 Package Insert

▲ Points to consider

- Use matched LAL and CSE reagents
- Plastic tubes are not recommended for making endotoxin dilutions
- Follow all suggested endotoxin vortexing times
- Use pyrogen-free accessories that have been qualified for endotoxin testing
- Equilibrate reagents to room temperature before use
- Do not vortex the LAL
- Avoid bubbles when plating reagents into the 96-well plate
- Avoid contaminating samples, standards, LRW and accessories
- Equipment should be installed, validated and maintained appropriately

Contact us

North America

Customer Service: +1800 638 8174 (toll free)

order.us@lonza.com

Scientific Support: +1800 521 0390 (toll free)

scientific.support@lonza.com

Europe

Customer Service: + 32 87 321 611

order.europe@lonza.com

Scientific Support: + 32 87 321 611 scientific.support.eu@lonza.com

International

Contact your local Lonza Distributor
Customer Service: +1301 898 7025
Fax: +1301 845 8291
scientific.support@lonza.com

Lonza Walkersville, Inc. – Walkersville, MD 21793

All trademarks belong to Lonza, registered in USA, EU or CH or to third party owners and used only for informational purposes. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights. The user bears the sole responsibility for determining the existence of any such third party rights, as well as obtaining any necessary licenses. For more details: www.lonza.com/legal. RT-DS013 01/22 ©2022 Lonza. All rights reserved.

www.lonza.com/kqcl www.lonza.com/turb