

**Kit SDS Cover Sheet****PyroCell™ MAT Cell Kit (M2016LC+M2016LS)**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

---

**Product Information**

Product name : PyroCell™ MAT Cell Kit (M2016LC+M2016LS)

Material number : 00249735

**Components**

MAT Cells

SDS attached

MAT Culture Medium Supplement

SDS attached

**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Cells**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

**SECTION 1. PRODUCT AND COMPANY IDENTIFICATION**

Product name : MAT Cells

**Manufacturer or supplier's details**Company : Capsugel Australia Pty Ltd.  
Suite 610, 12 Century Circuit,  
Norwest NSW 2153, AustraliaLonza Ltd  
Muenchensteinerstrasse 38  
CH-4002 Basel, Switzerland

Telephone : Tel +61 3 9550 0883

Telefax : Tel +61 3 9550 0890

E-mail address : sds@lonza.com

Emergency telephone number : +41 61 313 94 94 (24h)

**Recommended use of the chemical and restrictions on use**

Recommended use : For Research Use Only. Not for use in diagnostic procedures.

**SECTION 2. HAZARDS IDENTIFICATION****GHS Classification**

Not a hazardous substance or mixture.

**GHS label elements**

Not a hazardous substance or mixture.

**Other hazards which do not result in classification**

No information available.

**SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS**

Substance / Mixture : Mixture

**Hazardous components**

Chemical name	CAS-No.	Concentration (% w/w)
Substances with a workplace exposure limit		
Dimethyl sulfoxide	67-68-5	

**SECTION 4. FIRST AID MEASURES**

If inhaled : No special precautions required.

In case of skin contact : Wash with water and soap as a precaution.  
If skin irritation persists, call a physician.

**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Cells**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

---

In case of eye contact	:	Immediately flush eye(s) with plenty of water. If eye irritation persists, consult a specialist.
If swallowed	:	Immediately give large quantities of water to drink. Do not induce vomiting without medical advice. Never give anything by mouth to an unconscious person.
Most important symptoms and effects, both acute and delayed	:	No information available.
Notes to physician	:	No information available.

---

**SECTION 5. FIREFIGHTING MEASURES**

---

Suitable extinguishing media	:	Water spray Dry powder Foam
Hazardous combustion products	:	Carbon oxides (COx) Sulphur oxides (SOx)
Specific extinguishing methods	:	Use water spray to cool unopened containers.
Special protective equipment for firefighters	:	In the event of fire, wear self-contained breathing apparatus.

---

**SECTION 6. ACCIDENTAL RELEASE MEASURES**

---

Personal precautions, protective equipment and emergency procedures	:	Use personal protective equipment.
Environmental precautions	:	Prevent product from entering drains.
Methods and materials for containment and cleaning up	:	Contain spillage, and then collect with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations (see section 13).

---

**SECTION 7. HANDLING AND STORAGE**

---

Advice on protection against fire and explosion	:	Take precautionary measures against static discharges.
Advice on safe handling	:	No special precautions required.
Hygiene measures	:	Handle in accordance with good industrial hygiene and safety practice.

---

## SAFETY DATA SHEET

according to Work Health and Safety Regulations 2011

### MAT Cells

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

Conditions for safe storage

: No special storage conditions required.

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### Components with workplace control parameters

Components	CAS-No.	Value type (Form of exposure)	Control parameters / Permissible concentration	Basis
Dimethyl sulfoxide	67-68-5	TWA	250 ppm	WEEL

**Engineering measures** : Avoid splashes.

### Personal protective equipment

Respiratory protection : No personal respiratory protective equipment normally required.

Hand protection

Material : Nitrile rubber

Rate of permeability : > 480 min

Eye protection : Tightly fitting safety goggles

Skin and body protection : Choose body protection according to the amount and concentration of the dangerous substance at the work place.  
No special protective equipment required.

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : Aqueous solution

Colour : yellow-orange

Odour : no data available

Odour Threshold : no data available

pH : 6 - 8

Melting point/freezing point : no data available

Initial boiling point and boiling range : no data available

Flash point : no data available

Evaporation rate : no data available

Flammability (solid, gas) : no data available

**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Cells**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

---

Flammability (liquids)	:	no data available
Upper explosion limit / upper flammability limit	:	no data available
Lower explosion limit / Lower flammability limit	:	no data available
Vapour pressure	:	no data available
Relative vapour density	:	no data available
Relative density	:	no data available
Density	:	no data available
Water solubility	:	no data available
Partition coefficient: n-octanol/water	:	no data available
Auto-ignition temperature	:	no data available
Decomposition temperature	:	no data available
Viscosity, dynamic	:	no data available
Viscosity, kinematic	:	no data available

---

**SECTION 10. STABILITY AND REACTIVITY**

Reactivity	:	Stable under recommended storage conditions.
Chemical stability	:	Stable under normal conditions.
Possibility of hazardous reactions	:	Stable under normal conditions.
Conditions to avoid	:	no data available
Incompatible materials	:	Strong oxidizing agents
Hazardous decomposition products	:	No decomposition if stored normally.

---

**SECTION 11. TOXICOLOGICAL INFORMATION****Acute toxicity**

Acute oral toxicity	:	Remarks: no data available
Acute inhalation toxicity	:	Remarks: no data available
Acute dermal toxicity	:	Remarks: no data available

---

**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Cells**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

---

**Skin corrosion/irritation**

Remarks: no data available

**Serious eye damage/eye irritation**

Remarks: no data available

**Respiratory or skin sensitisation**

Remarks: no data available

**Germ cell mutagenicity**

Genotoxicity in vitro : Remarks: no data available

**Carcinogenicity**

Remarks: no data available

**Reproductive toxicity**

Effects on fertility : Remarks: no data available

**STOT - single exposure**

Remarks: no data available

**STOT - repeated exposure**

Remarks: no data available

**Further information**

Remarks: No data is available on the product itself.

---

**SECTION 12. ECOLOGICAL INFORMATION****Ecotoxicity**

Toxicity to fish : Remarks: no data available

**Persistence and degradability**

Biodegradability : Result: no data available

**Bioaccumulative potential**

Bioaccumulation : Remarks: no data available

**Mobility in soil**

Distribution among environmental compartments : Remarks: no data available

**Other adverse effects**

**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Cells**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

Additional ecological information : no data available

**SECTION 13. DISPOSAL CONSIDERATIONS****Disposal methods**

Waste from residues : Dispose of contents/container in accordance with local regulation.  
Contact waste disposal services.  
Solutions with high pH-value must be neutralized before discharge.

**SECTION 14. TRANSPORT INFORMATION****IATA** Not dangerous goods

**UN number** : Not applicable  
**Proper shipping name** : Not applicable  
**Transport hazard class** : Not applicable  
**Packing group** : Not applicable  
**Environmental hazards** : no

**IMDG** Not dangerous goods

**UN number** : Not applicable  
**Proper shipping name** : Not applicable  
**Transport hazard class** : Not applicable  
**Packing group** : Not applicable  
**Environmental hazards** : Marine pollutant: no

**ADR** Not dangerous goods

**UN number** : Not applicable  
**Proper shipping name** : Not applicable  
**Transport hazard class** : Not applicable  
**Packing group** : Not applicable  
**Environmental hazards** : no

## SAFETY DATA SHEET

according to Work Health and Safety Regulations 2011

### MAT Cells

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

<b>RID</b>	:	Not dangerous goods
<b>UN number</b>	:	Not applicable
<b>Proper shipping name</b>	:	Not applicable
<b>Transport hazard class</b>	:	Not applicable
<b>Packing group</b>	:	Not applicable
<b>Environmental hazards</b>	:	no
<b>Special precautions for user</b>	:	none
<b>Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code</b>	:	Not applicable

### SECTION 15. REGULATORY INFORMATION

**Safety, health and environmental regulations/legislation specific for the substance or mixture**

**Standard for the Uniform Scheduling of Medicines and Poisons**

Not applicable

### SECTION 16. OTHER INFORMATION

Revision Date : 30.06.2021

**Full text of other abbreviations**

WEEL : US. OARS. WEELs Workplace Environmental Exposure Level Guide, as amended

WEEL / TWA : Time weighted average

AiIC - Australian Inventory of Industrial Chemicals; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative)



**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Cells**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS - Workplace Hazardous Materials Information System

Date format : dd.mm.yyyy

**The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.**

AU / EN

## SAFETY DATA SHEET

according to Work Health and Safety Regulations 2011

### MAT Culture Medium Supplement

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

#### SECTION 1. PRODUCT AND COMPANY IDENTIFICATION

Product name : MAT Culture Medium Supplement

##### Manufacturer or supplier's details

Company : Capsugel Australia Pty Ltd.  
Suite 610, 12 Century Circuit,  
Norwest NSW 2153, Australia

Lonza Ltd  
Muenchensteinerstrasse 38  
CH-4002 Basel, Switzerland

Telephone : Tel +61 3 9550 0883

Telefax : Tel +61 3 9550 0890

E-mail address : sds@lonza.com

Emergency telephone number : +41 61 313 94 94 (24h)

##### Recommended use of the chemical and restrictions on use

Recommended use : For Research Use Only. Not for use in diagnostic procedures.

#### SECTION 2. HAZARDS IDENTIFICATION

##### GHS Classification

Not a hazardous substance or mixture.

##### GHS label elements

Not a hazardous substance or mixture.

##### Other hazards which do not result in classification

No information available.

#### SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

##### Hazardous components

Chemical name	CAS-No.	Concentration (% w/w)
Streptomycin sulphate	3810-74-0	>= 0,2 - < 0,3
Penicillin-G potassium salt	113-98-4	>= 0,1 - < 0,2

#### SECTION 4. FIRST AID MEASURES

If inhaled : Move to fresh air.  
Call a physician immediately.  
If breathing is irregular or stopped, administer artificial respiration.

**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Culture Medium Supplement**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

- 
- |   |   |   |
|---|---|---|
| In case of skin contact                                     | : | After contact with skin, wash immediately with plenty of soap and water.<br>Take off all contaminated clothing immediately.<br>Call a physician immediately.  |
| In case of eye contact                                      | : | Immediately flush eye(s) with plenty of water.<br>Call a physician immediately.   |
| If swallowed  | : | Immediately give plenty of water (if possible charcoal slurry).<br>Do not induce vomiting without medical advice.<br>Never give anything by mouth to an unconscious person.<br>Take victim immediately to hospital. |
| Most important symptoms and effects, both acute and delayed | : | No information available.   |
| Notes to physician  | : | No information available.   |
- 

**SECTION 5. FIREFIGHTING MEASURES**

- |   |   |  |
|---|---|--|
| Suitable extinguishing media                  | : | Water spray<br>Dry powder<br>Foam                              |
| Hazardous combustion products                 | : | No hazardous combustion products are known                     |
| Specific extinguishing methods                | : | Use water spray to cool unopened containers.                   |
| Special protective equipment for firefighters | : | In the event of fire, wear self-contained breathing apparatus. |
- 

**SECTION 6. ACCIDENTAL RELEASE MEASURES**

- |   |   |   |
|---|---|---|
| Personal precautions, protective equipment and emergency procedures | : | Use respirator when performing operations involving potential exposure to vapour of the product.  |
| Environmental precautions   | : | Do not flush into surface water or sanitary sewer system.   |
| Methods and materials for containment and cleaning up               | : | Contain spillage, and then collect with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations (see section 13). |
- 

**SECTION 7. HANDLING AND STORAGE**

- |   |   |  |
|---|---|--|
| Advice on protection against fire and explosion | : | Take precautionary measures against static discharges.         |
| Advice on safe handling                         | : | In case of insufficient ventilation, wear suitable respiratory |

## SAFETY DATA SHEET

according to Work Health and Safety Regulations 2011

### MAT Culture Medium Supplement

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

- equipment.  
 Avoid exposure - obtain special instructions before use.  
 Use only in area provided with appropriate exhaust ventilation.  
 Avoid contact with skin and eyes.
- Hygiene measures : DANGER!  
 Avoid contact with the skin and the eyes.  
 Avoid breathing dust or spray mist.  
 Wash hands before breaks and immediately after handling the product.
- Conditions for safe storage : No special storage conditions required.

## SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### Components with workplace control parameters

Contains no substances with occupational exposure limit values.

**Engineering measures** : Avoid splashes.

### Personal protective equipment

Respiratory protection : In the case of vapour formation use a respirator with an approved filter.

Hand protection

Material : Nitrile rubber

Rate of permeability : > 480 min

Eye protection : Tightly fitting safety goggles

Skin and body protection : Choose body protection according to the amount and concentration of the dangerous substance at the work place.  
 Complete suit protecting against chemicals

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : Aqueous solution

Colour : orange

Odour : no data available

Odour Threshold : no data available

pH : 6 - 8

Melting point/freezing point : no data available

Initial boiling point and boiling range : no data available

**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Culture Medium Supplement**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

---

Flash point	:	no data available
Evaporation rate	:	no data available
Flammability (solid, gas)	:	no data available
Flammability (liquids)	:	no data available
Upper explosion limit / upper flammability limit	:	no data available
Lower explosion limit / Lower flammability limit	:	no data available
Vapour pressure	:	no data available
Relative vapour density	:	no data available
Relative density	:	no data available
Density	:	no data available
Water solubility	:	no data available
Partition coefficient: n-octanol/water	:	no data available
Auto-ignition temperature	:	no data available
Decomposition temperature	:	no data available
Viscosity, dynamic	:	no data available
Viscosity, kinematic	:	no data available

---

**SECTION 10. STABILITY AND REACTIVITY**

Reactivity	:	Stable under recommended storage conditions.
Chemical stability	:	Stable under normal conditions.
Possibility of hazardous reactions	:	Stable under normal conditions.
Conditions to avoid	:	no data available
Incompatible materials	:	Strong oxidizing agents
Hazardous decomposition products	:	No decomposition if stored normally.

---

**SECTION 11. TOXICOLOGICAL INFORMATION**

**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Culture Medium Supplement**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

---

**Acute toxicity**

Acute oral toxicity : Remarks: no data available

Acute inhalation toxicity : Remarks: no data available

Acute dermal toxicity : Remarks: no data available

**Skin corrosion/irritation**

Remarks: no data available

**Serious eye damage/eye irritation**

Remarks: no data available

**Respiratory or skin sensitisation**

Remarks: no data available

**Germ cell mutagenicity**

Genotoxicity in vitro : Remarks: no data available

**Carcinogenicity**

Remarks: no data available

**Reproductive toxicity**

Effects on fertility : Remarks: no data available

**STOT - single exposure**

Remarks: no data available

**STOT - repeated exposure**

Remarks: no data available

**Further information**

Remarks: No data is available on the product itself.

---

**SECTION 12. ECOLOGICAL INFORMATION****Ecotoxicity**

Toxicity to fish : Remarks: no data available

**Persistence and degradability**

Biodegradability : Result: no data available

**Bioaccumulative potential**

Bioaccumulation : Remarks: no data available

**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Culture Medium Supplement**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

---

**Mobility in soil**

Distribution among environmental compartments : Remarks: no data available

**Other adverse effects**Additional ecological information : no data available

---

**SECTION 13. DISPOSAL CONSIDERATIONS****Disposal methods**Waste from residues : Dispose of contents/container in accordance with local regulation.  
Contact waste disposal services.  
Solutions with high pH-value must be neutralized before discharge.Contaminated packaging : Dispose of as unused product.

---

**SECTION 14. TRANSPORT INFORMATION****IATA** Not dangerous goods**UN number** : Not applicable  
**Proper shipping name** : Not applicable  
**Transport hazard class** : Not applicable  
**Packing group** : Not applicable  
**Environmental hazards** : no**IMDG** Not dangerous goods**UN number** : Not applicable  
**Proper shipping name** : Not applicable  
**Transport hazard class** : Not applicable  
**Packing group** : Not applicable  
**Environmental hazards** : Marine pollutant: no**ADR** :  
Not dangerous goods**UN number** : Not applicable  
**Proper shipping name** : Not applicable  
**Transport hazard class** : Not applicable  
**Packing group** : Not applicable  
**Environmental hazards** : no

## SAFETY DATA SHEET

according to Work Health and Safety Regulations 2011

### MAT Culture Medium Supplement

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

<b>RID</b>	:	Not dangerous goods
<b>UN number</b>	:	Not applicable
<b>Proper shipping name</b>	:	Not applicable
<b>Transport hazard class</b>	:	Not applicable
<b>Packing group</b>	:	Not applicable
<b>Environmental hazards</b>	:	no
<b>Special precautions for user</b>	:	none
<b>Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code</b>	:	Not applicable

#### SECTION 15. REGULATORY INFORMATION

**Safety, health and environmental regulations/legislation specific for the substance or mixture**

**Standard for the Uniform Scheduling of Medicines and Poisons**

Not applicable

#### SECTION 16. OTHER INFORMATION

Revision Date : 30.06.2021

##### Full text of other abbreviations

AIIIC - Australian Inventory of Industrial Chemicals; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical



**SAFETY DATA SHEET**

according to Work Health and Safety Regulations 2011

**MAT Culture Medium Supplement**

Version 2.0

Revision Date 30.06.2021

Print Date 24.07.2021

Substance Inventory; TDG - Transportation of Dangerous Goods; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS - Workplace Hazardous Materials Information System  
Date format : dd.mm.yyyy

**The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.**

AU / EN